

EdaMove 3

Electrodermal Activity and Physical Activity Sensor

Quick Start Guide

Welcome!

Thank you for choosing the EDA and physical activity sensor **EdaMove 3** from movisens. This sensor represents the latest technology for the mobile monitoring of electrodermal and physical activity in everyday life.

Scope of Application

The **EdaMove 3** captures the physical activity, EDA-parameters (i.e. SCL, SCR...) and other secondary parameters of a test subject.

It's been designed and tailored for research and scientific applications.

Software Installation

Before starting a measurement, you have to install the programs delivered on USB stick.

- **SensorManager:** Configure the sensor, start a measurement and download the data from the sensor.
- **UnisensViewer:** View and preprocess the stored data.
- **DataAnalyzer:** Analyze the stored sensor data (30 days free trial version).

Description of the Sensor

Featuring an integrated Bluetooth Smart Interface, the **EdaMove 3** expands researchers' options by allowing integration with the experience sampling platform **movisensXS**.

Comfortably worn on the wrist, the sensors electrodes attach to the palm via adhesive rings. A multi-colored LED provides visual feedback on the sensors operation, and a vibration alarm activates upon Bluetooth disconnection.




EdaMove 3

Electrodermal Activity and Physical Activity Sensor

Quick Start Guide

Handling



- **Charge the Sensor**
Connect the sensor to a USB port until the LED is blinking green. During charging, the status LED is blinking blue.
- **Start SensorManager**
Connect **EdaMove 3** to the PC, click on “Start Recording” and configure the measurement parameters.
- **Attach the Sensor**
Fill the electrodes with EDA-gel and place them with adhesive rings at the palms. Don't connect the sensor to a PC during measurement, this will stop the sensor.
- **Start SensorManager**
After the measurement, connect the **EdaMove 3** to the PC, click on „Save Data” and enter the required information.
- **Start UnisensViewer**
Or click on “Show data” directly after storing the data.
- **Start DataAnalyzer**
Select the desired analysis output to calculate secondary parameters or to generate reports.

